


## Photovoltaic Characterizations: Polarization and Mott Schottky plot

### I- Introduction

With the greenhouse effect and the increase of the price of the energy, new ways to produce energy are intensively developed. They include fuel cells, biofuel cells, nuclear, biomass, wind power, and photovoltaic. Among all these renewable energies, photovoltaic power seems very promising. Indeed, only 0.2% of the solar energy which touched the Earth surface will be sufficient to produce energy for all world [1,2]. A lot of work is performed in this field in order to improve the efficiency of such an energy supply.

The basic principle of a photovoltaic cell is presented in Fig. 1. It is constituted of two layers of *n*-type (electron acceptor) and *p*-type (electron donor) semiconductors.



**Fig. 1: Photovoltaic cell principle.**

In this note, several electrochemical investigations are performed in order to characterize the photovoltaic cell, such as I-V characterizations or electrochemical impedance spectroscopy (EIS).

### II- Experimental

Investigations are carried out with the SP-150 driven by EC-Lab<sup>®</sup> software. The size of the photovoltaic cell is 5.7 x 5.0 cm.


Measurements are performed under natural solar light during a sunny day.

**Note:** All settings and raw data files presented hereafter are available in the Data Sample folder of EC-Lab<sup>®</sup> Software as LP\_PV.mpr and Mott-Schottky\_PV.mpr.

### III- I-E characterization

Firstly, the photovoltaic cell is characterized by I-V characterization technique (Fig. 2).

Note that power is automatically ticked in the “Cell Characteristics” window if I-V characterization technique is selected in the “Photovoltaic/Fuel Cells” part.



**Fig. 2: “Parameters settings” window of I-V characterization technique.**

The resulting I-V characterization shows a typical *I* vs. *E* and *P* vs. *E* curves (Fig. 3). Several parameters can be drawn from this curve with the “Photovoltaic analysis” tool (Fig. 5):

- the Short Circuit Current ( $I_{sc}$ ), which corresponds to the maximum current when  $E = 0$  V,  $I_{sc} = 41$  mA,
- the Open Circuit Voltage ( $E_{oc}$ ), which is the potential when the current is equal to zero ampere,  $E_{oc} = 3.145$  V,
- the theoretical power ( $P_T$ ), which is defined by the following relationship  $P_T = I_{sc} \times E_{oc}$ ,  $P_T = 129$  mW,
- the maximum power,  $P_{MAX} = 91$  mW,

- the fill factor (FF), which is the ratio of  $P_{MAX}$  and  $P_T$ ,  $FF = 70.3\%$ ,
- the efficiency can also be calculated. If we assume that the solar power is equal to the means of the surfacic solar power ( $175 \text{ W/m}^2$ ), which is  $499 \text{ mW}$  for our photovoltaic cell surface area ( $28.5 \text{ cm}^2$ ). The efficiency of the solar cell is  $18\%$ .




Fig. 3: I-V characterization measurements.





Fig. 4: Photovoltaic analysis window.

## IV- Impedance characterization

### IV-1 Mott-Schottky

As photovoltaic cells are made of a semiconductor, the Mott-Schottky plot gives useful information [3,4]. This plot is available from Staircase Potentio Electrochemical Impedance Spectroscopy (SPEIS) investigation.

Settings are as follows (Fig. 5) and results are given Fig. 6:



(~ 1mn38s / scan)  
(dE ~ 0.100 V)

Fig. 5: SPEIS settings window.




Fig. 6: Bode plot of SPEIS investigation. Arrows show the shift of the modulus or phase shift during measurement.

The Mott-Schottky relationship involves the apparent capacitance measurement as a function of potential under depletion condition:

$$\frac{1}{C_{sc}^2} = \frac{2}{e\epsilon\epsilon_0 N} \left( E - E_{FB} - \frac{kT}{e} \right)$$

where  $C_{sc}$  is the capacitance of the space charge region,  $\epsilon$  is the dielectric constant of the semiconductor,  $\epsilon_0$  is the permittivity of free space,  $N$  is the donor density (electron donor concentration for an  $n$ -type semi-conductor or hole acceptor concentration for a  $p$ -type semi-conductor),  $E$  is the applied potential,  $E_{FB}$  is the flatband potential,  $k$  is the Boltzmann Constant, and  $T$  is the temperature.

The donor density can be calculated from the slope of the  $1/C^2$  vs.  $E_{we}$  curve, and the flatband potential can be determined by extrapolation to  $C = 0$ . The model required for the calculation is based on two assumptions:

a) Two capacitances must be considered: the one of the space charge region and the one of the double layer. Since capacitances are in series, the total capacitance is the sum of their reciprocals. As the space charge capacitance is much smaller than the double layer capacitance, the double layer capacitance contribution to the total capacitance is negligible. Therefore, the capacitance value calculated from this model is assumed to be the value of the space charge capacitance.

b) The equivalent circuit used in this model is a series combination of a resistor and a capacitance (the space charge capacitance). The capacitance is calculated from the imaginary component of the impedance ( $\text{Im}(Z)$ ) using the relationship:

$$\text{Im}(Z) = \frac{-1}{2\pi f C}$$

This model is adequate if the frequency is high enough (on the order of kHz). For the photovoltaic cell, the frequencies of interest are 200, 50, and 20 kHz (Fig. 6). For a frequency equal to 200 kHz,  $E_{FB}$  is 3.204 V, and the donor density is  $0.644 \cdot 10^{15} \text{ cm}^{-3}$ . The negative slope of the Mott-Schottky plot corresponds to a  $p$ -type conductivity of the photovoltaic cell.

The straight line in the Mott-Schottky plot is characteristic of the dopant with uniform distribution within the photovoltaic cell [4] (Fig. 7).





Fig. 7: Mott-Schottky plot at frequency of 200 (blue line), 50 (red line), and 20 kHz (green line).

#### IV-2 Equivalent circuit

The equivalent circuit of the photovoltaic cell can be figured by a  $R_1+R_2/C_2$  since one system is observed at around 30 kHz on the Bode plot (phase vs. frequency, in Fig. 6) [5]. Fits are performed at several applied potentials thanks to the "Z Fit" [6,7] tool. It is worth noting that fits on all cycles are

performed at once automatically by choosing "All" in the select option (Fig. 8).

Nyquist plots (Fig. 7) are typical for this kind of material [8].




Fig. 8: Nyquist plot.




Fig. 9: "Z Fit" tool.

The results of an EIS data fit are summarized in Fig. 10. As expected,  $R_1$ , which corresponds to ohmic resistance, is stable over the potential sweep.  $R_2$  decreases when the potential increases. It is important that  $R_2$  is also available from the  $I-E$  characterization. Indeed,

the inverse of the slope of the steady-state curve ( $\langle i \rangle$  vs.  $E_{we}$ , Fig. 3) is equal to this resistance. Some values are compared in Table 1. Both methods give similar results.

Table 1: Comparison between  $R_2$  determined by  $I-E$  and EIS characterization.

$E/V$	2	2.5	3
$R_2/\Omega$ from $I-E$ data	439	74	11
$R_2/\Omega$ from EIS data	368	86	9

$C_2$  is stable for a potential inferior to 2 V. Then decreases to OCV and finally increases.




Fig. 10: Zfit results with the following equivalent circuit  $R_1+R_2/C_2$ .

## V- Conclusion

This note exhibits the photovoltaic solar cells characterization by polarization and EIS techniques which allow the user to determine cell performance ( $I_{sc}$ ,  $E_{oc}$ ,  $P_{MAX}$ ,  $P_T$ , FF) and model (equivalent circuit).

Moreover, automatic fits on several EIS experiments are available in the "Z Fit" tool. More generally, this note shows the contribution of the electrochemistry in energy fields, which is currently a hot topic.

## References

- [1] Harnessing Solar Energy for the Production of Clean Fuel, T. Aartsma, E.-M. Aro, J. Barber, D. Bassani, T. Flüeler, H. de Groot, A. Holzwarth, O. Kruse, A. W. Rutherford, European Science Foundation, Science Policy Briefing, n°34, 2008.
- [2] Inorganic Photovoltaic Solar Cells: Silicon and Beyond, M. Tao, The Electrochemical Society Interface, Winter 2008, 30-35.
- [3] Diffusion Length Variation in Photovoltaic Cells with Bridgman-grown CuInSe<sub>2</sub> Substrates, C.H. Champness, Thin Solid Films, 2007, 515(15), 6200-6203.
- [4] On doubt about Mott-Schottky plot of organic planar heterojunction in photovoltaic cell, G. Jarosz, J. of Non-crystalline Solids, 2008, 354, 4338-4340.
- [5] Characterization of n-Si/RbAg<sub>4</sub>I<sub>5</sub> Interfaces by Photocurrent Measurements and Electrochemical Impedance Spectroscopy, N. Ohta, K. Takada, T. Sasaki, M. Watanabe, J. Electrochem. Soc., 2005, 152 (6), A1241-A1247.
- [6] Application Note #18, <http://www.bio-logic.info/potentiostat/notes.html>
- [7] Application Note #20, <http://www.bio-logic.info/potentiostat/notes.html>
- [8] Correlation between Photovoltaic Performance and Impedance Spectroscopy of Dye-Sensitized Solar Cells Based on Ionic Liquids, F. Fabregat-Santiago, J. Bisquert, E. Palomares, L. Otero, D. Kuang, S. Zakeerudin, M. Grätzel, . Phys. Chem. C 2007, 111, 6550-6560.